


# PediPlates<sup>®</sup>

I-Plate and O-Plate

*Indicated for angular  
deformity correction*

*Manufactured in  
stainless steel*

*Seven different sizes in  
2 plate configurations*



## Features and Benefits

- Stainless Steel Plates and Screws provide greater strength, resistance to breakage, as well as ease of removal
- Multiple plate sizes and configurations offer greater options intraoperatively
- All screws are self-tapping for easy insertion
- Utilizes 4.5 Stainless Steel Cannulated Screws with multiple screw length options
- Low-profile plates
- The I-Plate more adequately addresses larger and heavier patients by providing additional fixation options
- Dark laser etched line or center hole for guide wire on the plate allows for proper positioning



## PediPlates®



The O-Plate is offered in 4 sizes with 2 screw holes per plate.



The I-Plate is offered in 3 sizes with 4 screw holes per plate.



Before correction

During correction

The Adverse Effects and Contra-Indications listed in the IFU are important considerations. These Risks and general surgical Risks should be explained to the patient prior to surgery. †

† IFU, cleaning instructions, and surgical techniques may be obtained by calling OrthoPediatrics® Customer Service at 574-268-6379. Read and understand indications, warnings, and adverse effects explained in IFU's prior to use.

*OrthoPediatrics, Children Are Not Just Small Adults, PediPlates, PediLoc, PediFlex, PLEO, PediNail, Scwire, and the Pedi logo are trademarks of OrthoPediatrics Corp.*

*OrthoPediatrics, Children Are Not Just Small Adults, PediPlates, PediLoc, Scwire and the Pedi logo are registered trademarks in the United States. The Pedi logo is a registered Community Trade Mark.*



2850 Frontier Drive • Warsaw, IN 46582 • ph: 574.268.6379 or 877.268.6339 • fax: 574.268.6302 • [www.OrthoPediatrics.com](http://www.OrthoPediatrics.com)

OrthoPediatrics Corp. ©2013 Part# SS-1010-01-01 Rev A